



dermapose
Access


Guided fat harvest system to simplify office based procedures



www.hcp.dermapose.com

Discover What Is Possible
With Regenerative Fat

How does it work?



Dermapose Access seamlessly integrates to
Dermapose Refresh

Discover more on www.hcp.dermapose.com/access

ELEVATE SKIN UNDER VACUUM DIRECTLY INTO THE 90 ML ACCESS CAVITY. ARTICULATE INFILTRATION AND HARVESTING CANNULA WITH MM PRECISION.



Place Dermapose Access



Elevate skin using vacuum



Insert Needle Access
Through Dermis



Insert Refresh through
Needle Access point



EASY FAT HARVEST

Dermapose Access is a sterile, single use, closed system intended to assist the harvesting of autologous fat grafts from a patient for plastic and reconstructive surgery applications.

AN OVERALL SYSTEM

Once inserted, the harvesting syringe can be moved in a controlled manner making it possible to harvest fat directly into Dermapose Refresh.



Bimini Health Tech
420 Stevens Ave, Suite 220
Solana Beach, CA 92075

1-858-348-8050
info@BiminiHealthTech.com
www.BiminiHealthTech.com

Dermapose Access an innovative harvesting device



PRECISION

Controls both the depth and range of motion in harvesting procedure



AIDS INFILTRATION AND HARVEST

Pressure assisted tumescent infiltration



MINIMALLY INVASIVE

Access adipose tissue without a stitch or a cut



EASY TO USE

Allows your team to easily harvest up to 50cc of lipoaspirate just below the dermis



dermapose

You need more information?

info@dermapose.com

www.hcp.dermapose.com



BIMINI
HEALTHTECH

Bimini Health Tech
420 Stevens Ave, Suite 220
Solana Beach, CA 92075

1-858-348-8050

info@BiminiHealthTech.com

www.BiminiHealthTech.com

500-040-01-B